


V. Norman
January 21, 1963

DETERMINATION OF NO₂ IN CIGARETTE SMOKE

The determination of NO₂ in the gaseous phase of cigarette smoke by a colorimetric procedure involves the production of pink diazosulfanilic acid. The postulated color reaction is as follows:



Reagents:

1. Stock solution. 0.1 gm of N-(1-naphthyl)-ethylene diamine dihydrochloride in 100 ml of water.
2. Absorbing reagent. 5 gm of sulfanilic acid in almost 1 liter of H₂O containing 140 cc of glacial acetic acid. 20 ml stock solution added and diluted to 1 liter.
3. Standard solution. 0.0203 gm NaNO₂/liter. Equivalent of 10 ul/ of NO₂/1 ml.

Apparatus:

1. Smoking machine.
2. Smoke collection flask (55 ml).

88482240

Calibration Curve

0.5 - 5.0 ml portions of the standard solution are added to absorbing solution in 25 ml volumetric flasks and the solutions are diluted to 25 ml. with absorbing solution.

Color is allowed to develop for 15-20 minutes with occasional shaking and the optical density of the solutions is determined at 556 mμ against an absorbing solution blank.